

# Futar<sup>®</sup> D Occlusion

## It's The Hardest Thing We've Ever Done.

From the moment you try it, you'll see that Futar D performs like nothing else on the market. Upon application, Futar D provides a smooth, mousse-like consistency. In just seconds, it's hard enough to be trimmed with a bur or scalpel, and the excess can be snapped off in the mouth. Yet, Futar D is flexible enough not to shatter.



## Futar D from Roydent: A revolution in hardness.



Johnson City, TN 37604  
800-992-7767 | [www.roydent.com](http://www.roydent.com)

**ROYDENT** Dental Products  
The Quality Solution



Futar D Occlusion is applied straight from the cartridge, onto the occlusal surfaces

- 30 second work time and 90 second set time = 2 minutes
- Final set hardness of Shore D 43
- Mousse-like qualities
- Trims easily with a scalpel or bur



Futar D Occlusion bite registration can be checked in the mouth by breaking off the excess material

Item Number	Description
K11932	Pink (2-50ml cartridges & 12 mix tips)
K11932B	Pink-Bulk (12-50ml cartridges, tips not included)

YOUR ROYDENT DISTRIBUTOR:

R-FutarD 9/05



Johnson City, TN 37604  
800-992-7767 | [www.roydent.com](http://www.roydent.com)

**ROYDENT** Dental Products  
The Quality Solution