

Endo Storage Systems



Endo Storage Systems

- Mix or match in any combination to design the Endo Storage System that is right for your office
- All boxes are impervious to stains
- Resistant to chemicals and disinfecting solutions
- Can be sterilized up to 480°F

Item Number	Description
410STD	Semi Box 2000
450STD	Basic Box 2100
480STD	Mini Box 2100
401STD	Step Module (insert only)
402STD	Socket Module (insert only)
403STD	Mini Step Module (insert only)
404STD	Mini Socket Module (insert only)
481STD	Mini Box with Step Module
482STD	Mini Box with Socket Module

YOUR ROYDENT DISTRIBUTOR:



Johnson City, TN 37604
800-992-7767 | www.roydent.com

ROYDENT Dental Products
The Quality Solution

Interim Stand

Interim Stand

- For chairside storage of root canal instruments during endodontic treatment
- Resistant to chemical solutions
- Can be sterilized up to 480 F
- Approximate diameter = 2"
- Unique sponge material does not leave fibers on files



*Interim Stand is one stand, one refill sponge
Sponge refill is a pkg of 55 sponges*

Item Number	Description
495STD	Interim Stand
496STD	Refill Sponges-55

YOUR ROYDENT DISTRIBUTOR:



Johnson City, TN 37604
800-992-7767 | www.roydent.com

R-Bx/Stnd 2/04

