

Introducing 2Seal easymiX

Take the mess out of mixing

- ✓ Automatic 1:1 mix ratio
- ✓ Precise placement with 360° Flex Intra Oral Tip



2Seal just became 2easy:

Now available in a convenient auto-mix syringe with 360° Flex Intra Oral Tip

- Can be applied directly into the root canal system, or on a traditional mix pad
- Precise application even in canals with difficult access
- Easy to use

- Working time: 4 hours
- Setting time: 8 hours
- Radiopacity: 13.6 mm/mm Al
- Flow: 36 mm

- Film thickness: 26 μ m
- Shrinkage: 1.76%
- Solubility: 0.31%

 **ROYDENT**[™] roydent.com
D E N T A L P R O D U C T S 1.800.992.7767

2Seal easymiX

The quality of your sealer is critical for a biocompatible and long-term tight root canal filling. Epoxy amine-based sealers have been thoroughly tested in more than 90 scientific studies and well-proven in clinical applications.

Excellent biocompatibility


- Highly biocompatible
- Extruded material is easily resorbed and encapsulated by the surrounding tissue
- Antimicrobial

For all obturation methods

- Recommended for warm and cold obturation techniques
- Optimal viscosity, reaches all canal areas
- Excellent radiopacity

Dense and durable

- Outstanding flow characteristics
- Low shrinkage
- Low solubility
- High adhesion to gutta percha
- Long-term dimensional stability



2Seal is available in a standard 2-tube pack and now also as 2Seal easymiX in a convenient double-chamber mixing syringe for easy and fast application.



Product	Product Description	Item #
Starter Kit	<ul style="list-style-type: none"> • 2Seal easymiX Mixing Syringe, 15 g • 20 x 360° Flex Intra Oral Tips • 1 Mixing Pad • 1 Illustrated Tip Card • 1 Organizer 	102500R
Mixing Syringe Refill	• 2 x 2Seal easymiX Mixing Syringe, 15 g	102501R
Mixing Tips	• 40 x 360° Flex Intra Oral Tips	102502R
2Seal – Tubes	<ul style="list-style-type: none"> • 1 x 4 ml Paste A • 1 x 4 ml Paste B • 1 x Mixing Pad 	1020USR

Call your Roydent Representative to order today! 1.800.992.7767



608 Rolling Hills Drive • Johnson City, TN 37604 • roydent.com • 1.800.992.7767