



## Strong, yet flexible.

---

- Pre-applied silicone stops, color-coded by length: 21 mm — white, 25 mm — yellow
- Heat-tempered steel for stiffness
- Sharp edges twisted to the tip
- Noncutting tip follows canal curve

*Packaged 6 to a box*

# More sizes, more “control” with Roydent C-Files.

Size	21mm	25mm
6	06521006R	06525006R
8	06521008R	06525008R
10	06521010R	06525010R

Contact your local dealer



YOUR ROYDENT DISTRIBUTOR:



Johnson City, TN 37604  
800-992-7767 | [www.roydent.com](http://www.roydent.com)

R-C-File 2/04

**ROYDENT** Dental Products  
The Quality Solution