

Endodontic Forceps



R6661 – Advanced ME Forceps – 45°



R5606 – Micro Endo Forceps – 90°



R5716 – Post & Wedge Forceps

Endo Forceps

R5716 – Post & Wedge Forceps:

- Ideal for holding posts during placement
- Excellent for placing and removing wedges
- Concave and serrated jaws for firm grip
- Manufactured from German stainless steel

R5606 – Micro Endo Forceps – 90°:

- Serrated jaws for firm grip
- Ideal for removing silver points
- Manufactured from German stainless steel

R6661 – Advanced ME Forceps – 45°:

- Ideal for use during endodontic procedures in removing small objects from tight areas
- Advanced design offers a much shorter instrument with the beaks at a 45 degree angle for better access
- Manufactured from German stainless steel

YOUR ROYDENT DISTRIBUTOR:



Johnson City, TN 37604
800-992-7767 | www.roydent.com

R-EForceps 2/04

